

# Clifton Boiling Baths NE1B Series Boiling Baths

Clifton Range®



NE1B-14




Designed for applications that require continuous boiling. This range of baths can be regulated from gentle simmer to vigorous boil. When low water level is detected the red warning light on the front panel illuminates and heating is suspended. Once refilled heating is restored, and temperature control resumes automatically. With the heater directly immersed in water it provides even heat distribution across the tank base. The perforated tray ensures an even spread of heat throughout the chamber.

Comprising of a high quality corrosion resistant stainless steel tank with bright clean finish housed in a durable painted outer case. Painted surfaces feature anti-bacterial finish with good chemical resistance.

The analogue dial has a proportional 1-10 scale, and the energy regulator provides control from gentle simmer to vigorous boil.

	NE1B - 14	NE1B - 18W	NE1B - 28
<b>GENERAL</b>			
Temp. Range	Simmer to boiling at 100°C		
Temp. Control	Analogue dial, energy regulator		
Heater Power (220-240V)	2500W		
Water Level Alarm	Features low liquid level float device, heating automatically suspended		
Voltage	230V / 50-60Hz		
Drain (12mm Bore)	N/A	Drain Outlet	
<b>DIMENSIONS</b>			
Tank Capacity	14 Litre	18 Litre	28 Litre
External (mm)	362w x 332d x 290h		538w x 332d x 290h
Internal (mm)	325w x 300d x 150h	325w x 300d x 200h	500w x 300d x 200h



## Low liquid float device

A low liquid float device protects the heaters and a red warning light on the front panel alerts users to refill the bath.

# NE1B Series Boiling Baths

## Accessories



### Stainless Steel Gabled Lid with Insulated Handle

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	N/A	N/A	SL1-22SS



### Stainless Steel Flat One Piece Lid with Insulated Handle

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	LD-14		LD-22



### Stainless Steel Hinged Gabled Lid with Insulated Handle

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	SL1-14H/SS		N/A



### Clear Hinged Gable Lid

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	GL1-14		GL1-22



### Clear Lift Off Gable Lid

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	GL1-14AF		GL1-22AF



### Stainless Steel Lid with Concentric Rings

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	SLR1-14		SLR1-22
	SLR2-14		SLR3-22
	SLR3-14		

SLR1-4: 2 x 114mm holes

SLR1-14: 4 x 114mm holes

SLR1-22: 6 x 114mm holes

SLR1-8: 4 x 83mm holes

SLR2-14: 6 x 83mm holes

SLR3-22: 6 x 114mm holes

SLR2-8: 2 x 114mm holes

SLR3-14: 4 x 114mm holes



### Polypropylene Spheres (Pack of 200)

Bath	NE1B-14	NE1B-18W	NE1B-28
No. of Packs	2 Packs		4 Packs
Accessory Code	BP0368		

# NE1B Series Boiling Baths

## Accessories



### Raise and Lower Sterilising Tray

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	N/A	ST-1	N/A



### Pair of guide rail assemblies for ST-1 tray

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	N/A	SA01409	N/A



### Stainless Steel Test Tube Racks

Bath	NE1B-14	NE1B-18W	NE1B-28
No. of Racks	4		6
Accessory Code	6870 : 26 holes x 17mm DIA		
	6871 : 16 holes x 26mm DIA		
	6872 : 36 holes x 13mm DIA		
	6873 : 18 holes x 19mm DIA		
	6900 : 12 holes x 32mm DIA		
	6920 : 9 holes x 40mm DIA		



### Stainless Steel Flat n Fold™ Test Tube Racks

Bath	NE1B-14	NE1B-18W	NE1B-28
No. of Racks	4		6
Accessory Code	6770 : 26 holes x 17mm DIA		
	6771 : 16 holes x 26mm DIA		
	6772 : 36 holes x 13mm DIA		
	6773 : 18 holes x 19mm DIA		
	6800 : 14 holes x 32mm DIA		



### Drain Tap

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	N/A	BX0953	

## NE1B Series Boiling Baths Accessories



### Stainless Steel Wire Basket

Bath	NE1B-14	NE1B-18W	NE1B-28
Accessory Code	BSK-14		BSK-22

## Universal Accessories

TC-1	Thermometer Clip complete with Bent Stem Spirit Filled Thermometer
CAL	Calibration Record
BX0616	Draining Syphon
BKF-1	Bar Keepers Friend for cleaning water baths and removing dirt